Part no.	311-4122100-011
Product name	機器說明書/FORA 6 GTel/英法西./FCS/
Spec	L100*W150mm/高級畫刊紙/80頁/80P/黑/膠裝
Designer	Tsaiyi
Color	 K:100  K:20

FORA[®] 6 GTel

Cellular Multi-Functional Monitoring System
Système de surveillance cellulaire multifonctionnel
Sistema de monitorización multifuncional móvil

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Safety Information

Read the following *Safety Information* thoroughly before using the device.

- Use this device **ONLY** for the intended use described in this manual.
- Do **NOT** use accessories which are not specified by the manufacturer.
- Do **NOT** use the device if it is not working properly or damaged.
- This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only. Always consult your doctor to have the results interpreted.
- The blood glucose test strip can be used for the testing of newborns.
- The β -Ketone, total cholesterol and uric acid test strips shall **NOT** be used for the testing of newborns.
- Before using this device, read all instructions thoroughly and practice the test.
- Keep the device and testing supplies away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
- Do **NOT** use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the correct operation.
- Proper maintenance and storage are essential in ensuring the longevity of your device. If you are concerned about

the accuracy of the measurement, please contact the place of purchase or customer service representative for assistance.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your test results are lower or higher than usual, and you do not have symptoms of illness, first repeat the test. If you have symptoms or continue to get results higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood sample to perform a test. Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for

individuals experiencing a hyperglycaemic-hyperosmolar state, with or without ketosis. Please consult the healthcare professional before use.

- The measurement unit used for indicating the concentration of blood glucose can have mg/dL or mmol/L. The approximate calculation rule for conversion of mg/dL in mmol/L is:

mg/dL	Divided by 18	= mmol/L
mmol/L	Times 18	= mg/dL

For example:

1) $120 \text{ mg/dL} \div 18 = 6.6 \text{ mmol/L}$

2) $7.2 \text{ mmol/L} \times 18 = 129 \text{ mg/dL}$ approximately

Introduction

Intended Use

This system is intended for use outside the body (in vitro diagnostic use) by people with diabetes at home and by health-care professionals in clinical settings as an aid to monitor the effectiveness of diabetes control. It is intended to be used for the quantitative measurement of the blood glucose, β -ketone, total cholesterol and uric acid levels in whole blood. It should not be used for the diagnosis of diabetes.

Professionals may test with capillary and venous blood sample. Use only heparin for anticoagulation of whole blood.


Home use is limited to capillary blood from the finger tip and the approved sites.

Test Principle

With different types of test strips, your FORA 6 GTel Multi-functional Monitoring System allows you to measure the amount of blood glucose, β -Ketone, total cholesterol and uric acid levels in whole blood. The testing is based on the measurement of electrical current generated by the reaction of different substance with the reagent of the strip. The meter measures the current, calculates the blood glucose, β -Ketone, total cholesterol or uric acid, and displays the result.

The strength of the current produced by the reaction depends on the amount of the substance in the blood sample.

Product Overview

1. Test Strip Slot with Strip Indication Light
2. Power Button
3. Power Indicator
4. **Home** Button
5. Display Screen
6. Test Strip Ejector
7. Battery Compartment
8. Service Button
9. USB Port

Screen Display

- | | |
|---------------------|-------------------|
| 1. Measurement Mode | 6. Last Result |
| 2. Sync | 7. Date & Time |
| 3. Logbook | 8. Battery Symbol |
| 4. Settings | 9. RF Symbol |
| 5. Two-way Message | |

NOTE:

- Your meter display alternates between two kinds of results every 3 seconds, one being your last measurement result, and the other being your blood glucose 7-day average.
- If there are no records in memory, the meter displays "----".

Getting Started

Initial Setup


Step 1: Enter the Setting Mode

Power on the meter (with no test strip inserted). Tap Settings on the Main Page.

1. Setting the reminder alarm

Tap Alarm on the Settings page.

Your meter has four reminder alarms. Tap the time field to set the alarm time. Tap ON/OFF to enable/disable alarms.

2. Choosing a language


Tap Language on the Settings page. Select the appropriate language for use. Tap Back.

3. Setting the speaking volume


Tap Volume on the Settings page.

There are four (4) speaking volume levels to choose from. Select the desired speaking volume. Tap Back.


4. Setting the time

- Tap Time on the Settings page.
- Tap ▲ or ▼ to set the correct date. Press OK.
- Tap ▲ or ▼ to set the correct time. Press OK.

5. Setting the flight mode

Tap Flight mode on the Settings page to switch ON/OFF flight mode.

6. Setting the hi/low range for blood glucose measuring parameters

Tap Next on the bottom of Settings page to go to the next settings page. Tap BG Range.

Glucose:

- Three measuring modes are displayed on the screen: No Tag/Before Meal/After Meal. Tap each mode to set the hi/low range for each measurement.
- The hi/low range will affect the background color of your glucose test result.
- If your test result is higher than the normal range, the background will be red. If it is lower than the normal range, the background will be yellow. The background will be green when readings are within the set threshold.

Note:

- You may set it up to four reminder alarms.
- When the alarm goes off, the device will automatically turn on. Press the HOME button or the Power button the mute the alarm, or insert a test strip to start the measurement.
- If you do not mute the alarm, the device will beep for 2 minutes and then switch off.
- If the device is idle for 1 minute, it will enter power-saving mode. Press the HOME button to resume.
- If the device is idle for 3 minutes, it will automatically turn off.

Before Testing

Calibration

You must calibrate the meter every time you begin to use a new vial of β -ketone/ Total cholesterol/ Uric acid test strips by setting the meter with the correct code. Test results may be inaccurate if the code number displayed on the monitor does not match the number printed on the strip label or strip foil pack.

How to Code Your Meter (for β -ketone/ Total cholesterol/ Uric acid test)

1. Power on the meter and insert the code strip. Wait until the code number appears on the display.

Note:

Make sure the code numbers on display, code strip, and test strip vial or foil pack are the same. The code strip should be within the expiry date; otherwise, an error message may appear.

2. Remove the code strip and the meter will return to the home page. This tells you that the meter has finished coding and is ready for β -ketone/ Total cholesterol/ Uric acid testing.

Checking the Code Number

You need to make sure that the code number displayed on the meter matches the number on your test strip vial or foil

pack before you proceed. If it matches, you can proceed with your test. If the codes do not match, please stop testing and repeat the calibration procedure. If the problem persists, contact Customer Service for help.

NOTICE:



The codes used in this manual are examples only; your meter may display a different code.

WARNING:

- It is important to make sure that the LCD displayed code is the same as the code on your test strip vial or foil pack before testing. Failure to do so will get inaccurate results.
- If the LCD displayed code is not the same as the code on your test strip vial and the code number cannot be updated, please contact Customer Service for assistance.

Testing With Blood Sample

Test Strip Appearance

1. Blood Glucose
2. β -Ketone
3. Total Cholesterol
4. Uric Acid

1 2 3 4

Inserting a Test Strip

Insert the test strip into its slot.

Important!

The front side of test strip should face up when inserting the test strip. Test results might be wrong if the contact bar is not fully inserted into the test slot.


Important!

To reduce the chance of infection:

- Never share a lancet or a lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.

Preparing the Lancing Device


1. Remove the cap.

2. Insert a new lancet firmly into the white lancet holder cup.

3. Remove the protective disk on the lancet.

Hold the lancet firmly in place and twist off the protective disk.

4. Replace the cap until it snaps or clicks into place.

5. Rotate the dial to set the desired lancing depth.

6. Pull the cocking control out until the orange bar appears on the release button window.

Obtaining the Blood Sample

Please follow the suggestions below before obtaining a drop of blood:

- Wash and dry your hands before starting.
- Select the puncture site either on your fingertips or other body parts.
- Rub the puncture site for about 20 seconds before penetration.

❖ Blood from the fingertip

1. Press the lancing device tip firmly against the lower side of your fingertip.
2. Press the release button to prick your finger. A click indicates that the puncture is complete.

❖ **Blood from sites other than the fingertip (For Blood Glucose Test Strip Only)**

Alternative site testing (AST) is when individuals check their blood glucose levels using other areas of the body other than the fingertips. The FORA test strips allow AST to be performed on sites other than the fingertips. Please consult your health care professional before you begin AST.

Alternative site sample results may be different from fingertip sample results when glucose levels are changing rapidly (e.g., after a meal, after taking insulin, or during or after exercise).

We strongly recommend that you perform AST **ONLY** at the following times:


- During a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do **NOT** rely on test results at an alternative sampling site, but use samples taken from the fingertip, if any of the following applies:

- You think your blood sugar is low.
- You are not aware of symptoms when you become hypoglycemic.
- The results do not agree with the way you feel.
- After a meal.
- After exercise.
- During illness.
- During times of stress.

To obtain a blood sample from the alternative sites, please rub the puncture site for approximately 20 seconds.

1. Replace the lancing device cap with the clear cap.

2. Pull the cocking control out until the orange bar appears on the release button window.

Important!

- Do not use results from alternative site samples to calibrate continuous glucose monitoring systems (CGMS), or for insulin dose calculations.
- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.

- Avoid lancing the areas with obvious veins to avoid excessive bleeding.
- It is recommended to discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Performing A Test

1. Power on the meter.
(Do NOT insert a test strip before powering on the meter)
2. Tap Measure or insert the test strip to start a measurement.

3. Obtain a blood sample.

Use the pre-set lancing device to puncture the desired site. Wipe off the first appeared drop of blood with a clean cotton swab. Gently squeeze the punctured area to obtain another drop of blood. Be careful NOT to smear the blood sample.


The blood sample size of each test should be at least,

- blood glucose test: 0.5 μ L
- β -ketone test: 0.8 μ L
- total cholesterol test: 3.0 μ L
- uric acid test: 1.0 μ L





4. Apply the blood sample.

Move your finger to meet the absorbent hole of the test strip and the drop will be automatically be drawn onto the test strip. Remove your finger until the confirmation window is filled. The meter begins to count down. Do not remove your finger until you hear a beep sound.

5. Read your result.

The results of your test will appear after the meter counts down to 1. The results will be stored automatically in the meter memory.

6. Select your measuring mode.

- General Tests - any time of day without regard to time since the last meal.
- AC - no food intake for at least 8 hours.
- PC - 2 hours after a meal.

7. Data Transmission

The meter automatically transmits new data after each measurement.

Disposing Used Test Strip and Lancet

To remove the used test strip, simply push the **Test Strip Ejector** button upward to eject the used test strip. The device will automatically turn off after the test strip is removed. To remove the used lancet, remove the lancet from the lancing device after you have finished testing. Discard your used strip and lancet properly in a puncture-resistant container.

Important!

The used lancet and test strip may be biohazards. Please consult your health-care provider for proper disposal which complies with your local regulations.

Control Solution Testing

Our Control Solution contains a known amount of substance that reacts with test strips and is used to ensure your device and test strips are working together correctly.

Test strips, control solutions, or sterile lancets may not be included in the kit (please check the contents on your product box). They can be purchased separately.

Do a control solution test when:

- ✓ you suspect the device or test strips are not working properly.
- ✓ your test results are not consistent with how you feel, or if you think the results are not accurate.
- ✓ you have dropped or think you may have damaged the device.

To perform the control solution test, do the following:

1. Power on the meter.
(Do NOT insert a test strip before powering on the meter)
2. Tap Measure or insert the test strip to start a measurement.

3. Apply the control solution. Shake the control solution vial thoroughly before use. Squeeze out a drop and wipe it off, then squeeze another drop and place it on the tip of the vial cap. Hold the device to move the absorbent hole of the test strip to touch the drop. Once the confirmation window is filled completely, the device will begin counting down.

Note:

To avoid contaminating the control solution, do not directly apply the control solution onto a strip.

4. Read and compare the result. After counting down to 1, the test result of the control solution will appear on the display. Compare this result with the range printed on the test strip vial or individual foil pack and it should fall within this range. If the test result is out of range, read the instructions again and repeat the control solution test.

Note:

- The meter will detect the difference between control solution and blood samples automatically. It will automatically mark the result as a control solution test with "QC" display.
- Control solution test results are stored in the memory.
- The control solution range printed on the test strip vial or individual foil pack is for control solution use only. It is not the recommended range or reference values.
- Refer to the Maintenance section for important information about your control solutions.

Out-of-range results:


If you continue to get results that fall outside the range printed on the test strip vial, it means that the meter and strips might not be working properly. Please contact your local customer service or place of purchase for assistance.

Reviewing Test Results

Your device stores the 1000 most recent test results along with respective dates and times in its memory. To enter the meter memory, tap the Logbook icon on the Main Page.

To review all test results, do the following:

1. Tap Logbook icon on the Main Page
The first reading you see is the last test result along with date, time and the measurement mode.
2. Press ▲ or ▼ repeatedly to review other test results stored in the device.

To review the day-average test results, do the following:

1. Tap AVG tab
Your 7-, 14-, 21- day average results measured in No Tag mode will appear on the display. Tap ▲ or ▼ to review 28-, 60-, 90- day average results.
2. Tap Before Meal/After Meal to review 7-, 14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode.
3. Exit the meter memory
Press the Home button to return to the Main Page.

Note:

- If using the device for the first time, the “---” icon will appear when you recall the test results or review the average result. This indicates that there is no test result in the memory.
- Control solution results are NOT included in the day average.

Transferring Data

Data Transmission via Mobile Network


The meter could upload test results to a data system through a mobile network if your healthcare provider needs to collect the data remotely. Usually, the settings and the SIM card deployment should have already been placed in the meter when you receive it. The test results will be uploaded right after every measurement without a manual transmission. The data will be encrypted during transmission and only the service provider can decipher it. No personal identity will be revealed in the transmission payload.

Following the instructions of your care service provider should enable you to get connected with the service contact. When you use the meter for the very first time, it's highly recommended that you contact your service provider to confirm that the transmission to your account at the server side is correct and successful.

To learn more about the service procedure or transmission, please contact the customer service of your service provider. Data transmission is an optional function that depends on your service provider. You may perform measurements without this function.


1. To initiate transmission

- (a) The meter starts the transmission automatically after each measurement, sending the new data and previous data that has not been sent.
- (b) With the meter turned on, press Sync to initiate transmission for test results that have not been uploaded.

2. During transmission

When the arrow symbol shows up, the data transmission has started. The message “Upload OK” will appear when the transmission is finished.

3. After transmission


The uploaded results will change colors to indicate their statuses.

Note:

- While the meter is in transmission mode, it will be unable to perform a blood glucose test.
- Poor signal quality may be caused by busy service or no cellular reception nearby.
- If you receive a poor signal message, take the device to a spot with better reception and try again. Keep away from the microwave or other devices that may interfere with the signal.
- While the flight mode is turned on, the meter will not be able to transmit data.
- Do not remove or replace the SIM card arbitrarily. The data transmission might be interrupted due to improper placement of the card. The SIM card is dedicated to assigned servers and is not compatible with general-purpose phone sets. Please contact your service provider for further instructions on any transmission problems.

If Data Transmission Is Unsuccessful

If the Internet signal is unstable, your device may not be able to communicate with your destination account at that time. If the device cannot communicate with your account, the following the message will appear: "Poor signal"

Your results will be stored in the logbook and the upload will resume after the next measurement and when the signal is stable. You can also tab Sync on the screen to restart the upload when there is a reliable signal.

Data Security

The device will only transmit your measurement data in a secured and encrypted way to the data server that your service provider manages. No personally identifiable information will be covered and revealed. Please contact your service provider to learn more about their privacy policy regarding your measurement data.

When you use the meter for the very first time, it's highly recommended that you contact your service provider to confirm that the transmission to your account at the server side is correct and successful.

If you wish to cancel the data transmission function, please contact your service provider.

TWO-WAY MESSAGES

FORA 6 GTel is equipped with two-way message function to provide you with health tips and reminders. Types of messages include yes/no questions, multiple-choice questions and alert messages. You can click on the answer

icons to answer the questions.

After entering the message box, the screen displays the most recent message. You can click Previous or Next button to move to the previous or next message.

You can upload your answer onto the database and download new messages from the website. Press Home Button to return to the Main Page.

SERVICE BUTTON

The meter is equipped with a Service Button.

Pressing and holding the button for more than 3 seconds will send a Request Assistance to the care service provider.

For more information about the Service Button, please contact customer service or your healthcare professional.

Warning:

The Service Button is not intended for an emergency. Please DO NOT expect an instant response or immediate assistance by pressing this button. Your care service provider will define the service agreement with you.

Maintenance

Battery

Your meter comes with one 3.7V Li-ion rechargeable battery.

Low Battery Signal

The meter will display the message «Low battery» with a «Warning Symbol» to alert you when the meter power is getting low.

The power is not enough to do a test. Please recharge the battery immediately.

Recharging the battery

The power adapter can be used as the battery charger when you need to recharge the Li-ion battery. To recharge the 3.7V Li-ion battery, connect the AC adapter to the meter and a power source.

The power indicator will light up in red, which indicates the battery is undergoing a recharge. It should take approximately 5 hours to be fully charged. If the power indicator is blinking in red when connecting the AC adapter to a power source, this indicates the meter is in a pre-charging stage. If the power indicator starts blinking in green, this indicates the recharge is nearly finished. The power indicator will diminish when the recharge completes.

CAUTION:

The Li-ion battery must be recharged with the FORA 6 GTel meter by connecting the AC adapter to an electrical outlet. Only this 3.7V rechargeable Li-ion battery (battery model no. ICP553136ARM) can be recharged by FORA 6 GTel with the connection of the AC adapter. Any other kind of rechargeable battery is not allowed or it may damage the meter.

CAUTION
RISK OF EXPLOSION IF BATTERY IS REPLACED BY
AN INCORRECT TYPE.

DISPOSE OF USED BATTERIES ACCORDING TO THE
INSTRUCTIONS.

Note:

- FORA 6 GTel can use one Li-ion rechargeable battery **ONLY** (battery model no. ICP553136ARM).
- Wrong types of batteries may cause an explosion.
- Do not use other types of batteries, mixing old and new or batteries of different brands. Use only new batteries of the specified size and type.
- Replacing the batteries does not affect the test results stored in the memory.
- Keep the battery away from small children. If swallowed, promptly seek medical assistance.
- Batteries might leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period.
- Properly dispose of the batteries according to your local environmental regulations.

Caring for Your Device

- To clean the exterior of the device, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft dry cloth. Do **NOT** rinse with water.
- Do **NOT** use organic solvents to clean the device.

Device Storage

- Storage condition: -20°C to 60°C (-4°F to 140°F), below 95% relative humidity.
- Always store or transport the device in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

Meter Disposal

The used meter should be treated as contaminated and may carry a risk of infection during measurement. The batteries in this used meter should be removed and the meter should be disposed in accordance with local regulations.

Caring for Your Test Strips

- Storage conditions: 2°C to 30°C (35.6°F to 86°F) and below 85% relative humidity for all test strips (blood glucose, β -ketone, total cholesterol, uric acid). Do **NOT** freeze.
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool and dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.

- Do not use test strips beyond the expiry date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for assistance.

For further information, please refer to the test strip package insert.

Important Control Solution Information

- Use only our control solutions with your device.
- Do not use the control solution beyond the expiry date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test be done at room temperature 20°C to 25°C (68°F to 77°F). Make sure your control solution, device, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 2°C to 30°C (35.6°F to 86°F). Do **NOT** freeze.

DETAILED INFORMATION

The meter provides you with plasma equivalent results.

Desirable ranges:

Normal plasma glucose range for people without diabetes * ₁	Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
	2 hours after meals	< 140 mg/dL (7.8 mmol/L)
Total cholesterol * ₂		< 200 mg/dL (5.17 mmol/L)
Uric Acid * ₃	Male	3.5 ~ 7.2 mg/dL (0.208 ~ 0.428 mmol/L)
	Female	2.6 ~ 6 mg/dL (0.155 ~ 0.357 mmol/L)

- The β -Ketone test measures Beta-Hydroxybutyrate (β -OHB), the most important of the three β -Ketone bodies in the blood. Normally, levels of β -OHB are expected to be less than 0.6 mmol/L.
- β -OHB levels may increase if a person fasts, exercises vigorously or has diabetes and becomes ill. If your β -Ketone result is 0.0 mmol/L, repeat the β -Ketone test with new test strips. If the same message appears again or the result does not reflect how you feel, contact your healthcare professional.

Follow your healthcare professional's advice before you make any changes to your diabetes medication programme. If your β -Ketone result is between 0.6 and 1.5 mmol/L, this may indicate development of a problem that could require medical assistance. Follow your healthcare professional's

instructions. If your β -Ketone result is higher than 1.5 mmol/L, contact your healthcare professional promptly for advice and assistance. You may be at risk of developing diabetic ketoacidosis (DKA).

Please consult your doctor to determine a target range that works best for you.

- *1. American Diabetes Association (2014). Clinical Practice Recommendations. Diabetes Care, 37 (Supplement 1): S16.
- *2. Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III) final report(2002).Circulation106: 3143–3421
- *3. National Kidney Foundation (2014)
- *4. NIH - National Heart, Lung, and Blood Institute (NHLBI)
<https://www.nhlbi.nih.gov/health-topics/blood-tests>

Symbol Information

Symbol	Referent	Symbol	Referent
	For in vitro diagnostic use		Do not reuse
	Consult instructions for use		Temperature limitation
	Use by		CE Mark

	Batch code		Manufacturer
	Keep away from sunlight		Dispose of the packaging properly after use
	Serial number		Caution
	Keep Dry		Sterilized using irradiation
	Model number		Do not use if package is damaged
	Dispose of or recycle the electrical wastes according to local regulations		Authorised representative in the European Union
1.5V 	1.5 Volts DC		Unique Device Identification
	Quantity		Battery
	Humidity limitation		Appareil médical
	Importer		Distributor

Troubleshooting

If you follow the recommended steps but the problem persists, or error messages other than the ones below appear, please call your local customer service.

Do not attempt to repair the device yourself and never try to disassemble the device under any circumstances.

Result Readings

Blood Glucose Test:

Message	What it Means
 Low	< 10 mg/dL (0.5 mmol/L)
 High	> 600 mg/dL (33.3mmol/L)

β-Ketone Test:

Message	What it Means
 Low	< 0.1 mmol/L
 High	> 8.0 mmol/L

Total Cholesterol Test:

Message	What it Means
 Low	< 100 mg/dL (2.65 mmol/L)
 High	> 400 mg/dL (10.4 mmol/L)

Uric Acid Test:

Message	What it Means
 Low	< 3 mg/dL (0.179 mmol/L)
 High	> 20 mg/dL (1.190 mmol/L)

Error Message

Error Message	What to do
No strip inserted	Repeat the test with a new test strip. If the meter still does not work, please contact the customer service for assistance.
Used strip	
Strip incorrect 1-2	
Measuring cancelled	Review the instructions and repeat the test with a new test strip. Please make sure you use the correct technique and dose with enough blood..
Not enough blood	

No blood sample detected	Check if the absorbent hole is filled with blood or control solution. If the problem persists, please contact the customer service for assistance.
Temperature too high	The operating temperature range is from 8°C to 45°C (46.4°F to 113°F). Repeat the test after the device and test strip have reached the above-mentioned operating temperature.
Temperature too low	
BGM no response	Restart the meter. If the meter still does not work, please contact the customer service for assistance.
IC1 incorrect	
Operation Error	
SIM card not detected	Check if the SIM card is correctly inserted.
SIM card unaccepted 1-2	The SIM card cannot be used with the meter. Please contact your local customer service for assistance.
Unknown code card	Check if you have the correct code card. Repeat the coding procedure with the correct code card. If the problem persists, contact your local customer service for assistance.
Code card incorrect	
Account not exist	Register your meter online.

Server setup error	Please contact the customer service for assistance.
Poor signal 1-4	Press Sync button to upload data again when the signal is good. Or wait for the meter to automatically upload data again after the next measurement.
Invalid operation during Bat. charging	Remove the charging adapter and retry.
Low battery	
BGPRs Bat. not enough	
Unknown SW Err	Please contact the customer service for assistance.
Unknown Error	
Please turn off flight mode	Turn off Flight mode on the Settings page.

Blood Glucose Measurement

Symptom	Cause	What To Do
The device does not display a message after inserting a test strip.	Batteries exhausted.	Replace the battery immediately and reset date and time on the meter setting.
	Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
	Defective device or test strips.	Please contact customer service.
The test does not start after applying the sample.	Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
	Defective test strip.	Repeat the test with a new test strip.
	Sample applied after the device is automatically turned off.	Repeat the test with a new test strip. Apply sample only when the meter displays "Take a blood sample."
	Defective device.	Please contact customer service.

The control solution testing result is out of range.	Error in performing the test.	Read instructions thoroughly and repeat the test again.
	Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
	Expired or contaminated control solution.	Check the expiration date of the control solution.
	Control solution that is too warm or too cold.	Control solution, device, and test strips should be at room temperature (20°C to 25°C / 68°F to 77°F) before testing.
	Defective test strip.	Repeat the test with a new test strip.
	Device malfunction.	Please contact customer service.

Specifications

Memory	1000 measurement results with respective date and time
Dimensions	98 (L) x 55 (W) x 15.5 (H) mm
Power Source	One 3.7V Li-Ion rechargeable battery
Weight	63.2g (without battery)
External output	3G/4G Mobile network
Features	Auto electrode insertion detection Auto sample loading detection Auto reaction time count-down Auto switch-off after 3 minutes without action Temperature warning
Operating Condition	8°C to 45°C (46.4°F to 113°F), below 85% R.H. (non-condensing)
Storage/ Transportation Condition	-20°C to 60°C (-4°F to 140°F), below 95% R.H
Measurement Units	Blood Glucose Test: mg/dL or mmol/L β-ketone Test: mmol/L Total cholesterol Test: mg/dL Uric acid Test: mg/dL

Measurement Range	Blood Glucose Test: 10 ~ 600 mg/dL (0.55 ~ 33.3 mmol/L) β-ketone Test: 0.1 ~ 8.0 mmol/L Total cholesterol Test: 100 ~ 400 mg/dL (2.6 ~ 10.4 mmol/L) Uric acid Test: 3 ~ 20 mg/dL (0.179 ~ 1.190 mmol/L)
Hematocrit range	Blood Glucose Test: 0 ~ 70% β-ketone Test: 10 ~ 70% Uric Acid Test: 20 ~ 60% Total Cholesterol Test: 20 ~ 60%
Test Sample	Blood Glucose Test: Capillary/ Venous/ Neonatal/ Arterial β-ketone Test: Capillary/ Venous Uric Acid Test: Capillary Total Cholesterol Test: Capillary
Test Result	Measurements are reported as plasma equivalents

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-, IEC/EN 61326-2-6, EN 301 489-17, EN 300 328.

WARRANTY TERMS AND CONDITIONS

With respect to disposable products, ForaCare Suisse warrants to the original purchaser that, at time of delivery, each standard product manufactured by ForaCare Suisse shall be free from defects in material and workmanship and, when used for the purposes and indications described on the labeling, is fit for the purposes and indications described on the labeling. All warranties for a product shall expire as of the product expiration date, or if none, after two (2) year from the original date of purchase, as long as it has not been modified, altered, or misused. ForaCare Suisse warranty hereunder shall not apply if:

(i) a product is not used in accordance with its instructions or if it is used for a purpose not indicated on the labeling; (ii) any repairs, alterations or other work has been performed by the buyer or others on such item, other than work performed with ForaCare Suisse's authorisation and according to its approved procedures; or (iii) the alleged defect is a result of abuse, misuse, improper maintenance, accident or the negligence of any party other than ForaCare Suisse. The warranty set forth herein is conditioned upon proper storage, installation, use and maintenance in accordance with applicable written recommendations from ForaCare Suisse.

The warranty furnished hereunder does not extend to damaged items purchased hereunder resulting in whole or in part from the use of components, accessories, parts or supplies not furnished by ForaCare Suisse.

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- *1. American Diabetes Association (2014). Clinical Practice Recommendations. Diabetes Care, 37 (Supplement 1): S16.
- *2. Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III) final report(2002).Circulation106: 3143–3421
- *3. National Kidney Foundation (2014)
- *4. NIH - National Heart, Lung, and Blood Institute (NHLBI) <https://www.nhlbi.nih.gov/health-topics/blood-tests>

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ForaCare Suisse AG
MedNet EC-REP GmbH